




Sistema de Anclaje para Sobredentaduras
sobre implantes dentales



Distribuidor: BINTAL UNIVERSA, S.L.
PLANTIO 22 – POZUELO DE ALARCON – MADRID
91 352 22 99
bintal@bintaluniversal.com

TIPO IMPLANTE




Los pilares KERATOR se pueden utilizar para la adaptación a todos los sistemas de implantes. (Branemark, 3i, Astra, ITI, Osstem, Zimmer...)

VENTAJAS

1. Fácil inserción del pilar KERATOR por incluir un novedoso posicionador.
2. Altura vertical más baja del mercado.
3. Doble retención externa e interna que garantiza un rendimiento más duradero.
4. Mejor retención resiliente con mayor movilidad, facilidad de uso para el paciente.
5. Los pilares KERATOR se pueden utilizar con implantes divergentes, utilizando el exclusivo **Pilar KERATOR Angulado** o también el corrector de divergencia de 40° (incluido en el Set de Prótesis).

Pilar Angulado KERATOR



Utilizando nuestro propio diseño patentado del Pilar Angulado KERATOR, podemos cambiar la trayectoria o divergencia de los implantes. Maximizamos la satisfacción de los pacientes y odontólogos.

Cápsula de metal



Transportador



Con el novedoso posicionador, realizaremos una fácil inserción del pilar manualmente.

Componentes KERATOR



Manual KERATOR

Pilares Rectos y Angulados



Kit de herramientas
KERATOR



Extracción del pilar de
cicatrización



1- Pilar Angulado
2- Conector Angulado



Atornillar conector
angulado en el implante
con su transportador



Roscar conector



Ajustar el conector con
destornillador hex.
1.25mm con torque




Colocar pilar recto
KERATOR



Atornillar pilar recto
KERATOR



El conector angulado se
ajusta con la herramienta



Ajustar el pilar angulado
utilizando un martillo



Modelo terminado



Envase set de prótesis



Colocar el aro espaciador
en el pilar KERATOR



Colocar la cápsula de
metal con el retenedor
negro



Preparar unos huecos en
la prótesis



Resinar la prótesis



Colocar la prótesis en la
posición deseada



Extraer la prótesis y el
espaciador



Sustituir el retenedor
negro



Colocar la retención
definitiva



Prótesis terminada



Colocar prótesis

HEXED(External Type)		KERATOR IMPLANT ATTACHMENT PRODUCT CODES							
		HEIGHT OF ABUTMENT TISSUE CUFF ABOVE THE IMPLANT							
IMPLANT DESCRIPTION		0 mm	1.0 mm	2.0 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm	Angled
NOBEL BIOCARE	BRANEMARK 3.3 (NP)(HEX)								
OSSTEM	US _ 3.5, SS _3.5								
DIO	FSN 3.3		B331	B332	B333	B334	B335	B336	B33AN
MEGAGEN	EXFEEL EXTERNAL 3.3								
NEOBIOTECH	SINUS QUICK EB 3.5								
BIO-HORIZON	MAESTRO 3.5 (yellow) (HEX)		BHM351	BHM352	BHM353	BHM354	BHM355	BHM356	BHM35AN
	MAESTRO 4.0 (green) (HEX)		BHM401	BHM402	BHM403	BHM404	BHM405	BHM406	BHM4AN
	MAESTRO 5.0 (blue) (HEX)		BHM501	BHM502	BHM503	BHM504	BHM505	BHM506	BHM5AN
BIO-LOK	3.45 (green) (HEX) & Internal Connection		BL401	BL402	BL403	BL404	BL405	BL406	BL4AN
	5.0 (red) (HEX) & Internal Connection		BL501	BL502	BL503	BL504	BL505	BL506	BL5AN
NOBEL BIOCARE	HEX-LOCK (HL) 3.25 (HEX)		SH3251 (1.5mm)	SH3252 (2.5mm)	SH3253 (3.5mm)	SH3254 (4.5mm)	SH3255 (5.5mm)	SH3256 (6.5mm)	SH325AN
	REPLACE 3.5 (HEX)		SR351	SR352	SR353	SR354	SR355	SR356	SR35AN
ACE SURGICAL	4.0, 5.0 (4.1) (HEX)								
ATTACHMENTS INTERNATIONAL	IMZ 3.3 (HEX)								
	IMZ 4.0 (HEX)								
BIO-LOK	4.0 (yellow) (HEX) & Internal Connection								
INNOVA	4.1 ENDOPORE (.7mm HEX)								
	3.25 ENTEGRA (.7mm HEX)								
	4.1 PLATFORM 03B								
	4.0 ENTEGRA (.7mm HEX)								
IMPLANT INNOVATIONS (3I)	3.75 MINIPLANT (HEX)								
	4.0 CYLINDRICAL (HEX)								
IMTEC CORPORATION	SELF TAPPING 3.75 (HEX)								
	UNIVERSAL 3.4, 4.0 (HEX)								
INTRA-LOCK SYSTEM INTERNATIONAL	3.75, 4.0 HEX THREADED								
LIFECORE BIOMEDICAL	RESTORE 3.75, 4.0 (HEX)								
	SUSTAIN 4.2 (HEX)								
MIS MEDICAL IMPLANT SYSTEM	EXTERNAL HEXAGON (STANDARD)		B401	B402	B403	B404	B405	B406	B4AN
NOBEL BIOCARE	BRANEMARK 3.75, 4.0, 5.0 (RP)(HEX)								
	HEX-LOCK (HL) 3.8, 4.5 (HEX)								
OSTEO-IMPLANT CORPORATION	3.25 SD, 3.75 ST, 4.0 DT								
STERNGOLD IMPLAMED	HEX CYLINDER 3.3, 4.0								
	HEX SCREW 3.75, 4.0, 5.0								
ZIMMER - CALCITEK	SWEDE-VENT 4.0 (HEX)								
	SWEDE-VENT TL (HEX)								
	TAPERLOCK (HEX)								
	THREADLOCK 4.0 (HEX)								
OSSTEM	US _ 3.75, 4.1								
DIO	FSN 4.0								
MEGAGEN	EXFEEL EXTERNAL 4.0								
COWELL MEDI	ATLAS EXTERAL 4.0								
NEOBIOTECH	SINUS QUICK EB 4.0								
INNOVA	5.0 ENDOPORE (.9mm HEX)		IN501	IN502	IN503	IN504	IN505	IN506	IN5AN
ZIMMER - CALCITEK	SPLINE 3.25		CS3251	CS3252	CS3253	CS3254	CS3255	CS3256	CS325AN
	SPLINE 3.75, 4.0		CS401	CS402	CS403	CS404	CS405	CS406	CS4AN
	SPLINE 5.0		CS501 (.86 mm)	CS502	CS503	CS504	CS505	CS506	CS5AN
KYOCERA POI	KYOCERA POI	KP42035 (.35 mm)	KP421	KP422	KP423	KP424	KP425	KP426	KP42AN
SARGON	3.8 (HEX)		SI382 (1.8 mm)		SI383	SI384	SI385	SI386	SI38AN

INTRA-LOCK SYSTEM INTERNATIONAL	4.75 HEX THREADED		LC501 (1.42mm)	LC502	LC503	LC504	LC505	LC506	LC50AN
	5.0 THREADED CONIC								
LIFECORE BIOMEDICAL	LIFECORE 5.0 (HEX)								
LIFECORE BIOMEDICAL	RESTORE 3.3 (SD) (HEX)		LR331	LR332	LR333	LR334	LR335	LR336	LR33AN
BICON , ALLODEN	2.0 POST			BC202	BC203	BC204	BC205	BC206	BC2AN
	3.0 POST			BC302	BC303	BC304	BC305	BC306	BC3AN
NTRA-LOCK SYSTEM INTERNATIONAL	4.3 THREADED CONIC		SR431	SR432	SR433	SR434	SR435	SR436	SR43AN
NOBEL BIOCARE	REPLACE 4.3 (HEX)								
IMPLANT INNOVATIONS (3i)	3.25 MICRO MINI (HEX)								
	3i 5.0 (HEX)								
MIS MEDICAL IMPLANT SYSTEM	EXTERNAL HEXAGON (WIDE)		I501 (.6 mm)	I502	I503	I504	I505	I506	I5AN
DIO	FST 5.0								
NEOBIOTECH	SINUS QUICK EB 5.0 (3i Type)								
NEOBIOTECH	EXTERNAL WIDE 5.0								
NOBEL BIOCARE	BRANEMARK 5.0 (WP) (HEX)								
OSSTEM	US _ 5.0, 5.5		B501	B502	B503	B504	B505	B506	B5AN
DIO	FSN 5.0, 5.5								
MEGAGEN	EXFEEL EXTERNAL 5.0, 5.5								
NEOBIOTECH	SINUS QUICK EB 5.0 6.0 (BRANEMARK Type)								
NOBEL BIOCARE	HEX-LOCK (HL) 5.0 (HEX)								
NOBEL BIOCARE	REPLACE 5.0 (HEX)		SR501 (1.37 mm)	SR502	SR503	SR504	SR505	SR506	SR5AN
INNOVA	3.5, 4.0 ENDOPORE (.9mm HEX)		IN351	IN352	IN353	IN354	IN355	IN356	IN35AN
	3.25 ENTEGRA (.9mm HEX) 3.5 PLATFORM 03M								
IMTEC CORPORATION	MINI 3.3 (HEX)		IC331 (.91 mm)	IC332	IC333	IC334	IC335	IC336	IC33AN
NON-HEXED(Internal Type)	KERATOR IMPLANT ATTACHMENT PRODUCT CODES								
	HEIGHT OF ABUTMENT TISSUE CUFF ABOVE THE IMPLANT								
IMPLANT DESCRIPTION		0 mm	1.0 mm	2.0 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm	Angled
ATTACHMENTS INTERNATIONAL	IMZ 3.3 (NON-HEX)	IM330	IM331	IM332	IM333	IM334	IM335	IM336	IM33AN
	IMZ 4.0 (NON-HEX)	IM400	IM401	IM402	IM403	IM404	IM405	IM406	IM4AN
LIFECORE BIOMEDICAL	RENOVA 3.75	RN3750	RN3751	RN3752	RN3753	RN3754	RN3755	RN3756	RN375A
	RENOVA 4.5, 4.75	RN4750	RN4751	RN4752	RN4753	RN4754	RN4755	RN4756	RN475A
BIO-HORIZON	PRODIGY 5.7 (green) (Internal)		BHP571	BHP572	BHP573	BHP574	BHP575	BHP576	BHP57AN
ZIMMER - CALCITEK	INTEGRAL 3.25 SD	CI3250	CI3251	CI3252	CI3253	CI3254	CI3255	CI3256	CI325AN
	OMNIOLOC 3.25								
IMPLANT INNOVATIONS (3i)	TG		TG401	TG402	TG403	TG404	TG405	TG406	TG40AN
MIS MEDICAL IMPLANT SYSTEM	INTERNAL HEXAGON (STANDARD)								
ASTRA TECH	3,5		AS3508 (.81mm)	AS352	AS353	AS354	AS355	AS356	AS35AN
M	REGULAR , Mini								
Oneplant	INTERNAL 4.3 5.3	WA400	WA401	WA402	WA403	WA404	WA405	WA406	WA40AN
LIFECORE BIOMEDICAL	SUSTAIN 4.0, 4.7	CC400	CC401	CC402	CC403	CC404	CC405	CC406	CC4AN
NOBEL BIOCARE	STERIOSS 3.8								
ZIMMER - CALCITEK	INTEGRAL 4.0								
NOBEL BIOCARE	STERIOSS 3.25	SO3250	SO3251	SO3252	SO3253	SO3254	SO3255	SO3256	SO325AN
OSSTEM	GS _ , GS III 3.5		G351	G352	G353	G354	G355	G356	G3AN
OSSTEM	GS _ , GS III 4.0, 4.5, 5.0		G401	G402	G403	G404	G405	G406	G4AN
New ASTRA (Aqua)	1.6 Thread		AT401	AT402	AT403	AT404	AT405	AT406	AT4AN
MEGAGEN	EZ PLUS 3.5								
COWELL MEDI	Submerged								
ASTRA TECH	4.0, ST		AS401 (.5 mm)	AS402	AS403	AS404	AS405	AS406	AS4AN
DENTIUM	IMPLANTUM 3.4, 3.8, 4.3, 4.8 /SUPERLINE								
MEGAGEN	EZ PLUS 4.0 4.5 5.0								
	MAGAFIX 3.4 3.8 4.3 4.8 (HEX Type)								
DENTIS	DENTIS SUBMERGED 3.7 4.3 4.8								
DIO	SM REGULAR WIDE 4.5/5.3(PLATFROM 4.5mm)	D451	D452	D453	D454	D455	D456	D45AN	

DIO	SM NARROW 3.8		D381	D382	D383	D384	D385	D386	D38AN
BIO-HORIZON	PRODIGY 3.5 (yellow) (Internal)								
ZIMMER - PARAGON (3.5mm)	BIO-VENT 3.5	PA350	PA351	PA352	PA353	PA354	PA355	PA356	PA35AN
	CORE-VENT 4.3								
	MICRO-VENT 3.25								
	SCREW-VENT 3.3, 3.7								
ZIMMER - PARAGON	ADVENT 4.5		AI451	AI452	AI453	AI454	AI455	AI456	AI45AN
ZIMMER - CALCITEK	OMNILOC 4.0	CO400	CO401	CO402	CO403	CO404			CO4AN
BIO-HORIZON	PRODIGY 4.0 (green) (Internal)								
ZIMMER - PARAGON (4.5mm)	BIO-VENT 4.5	PA450	PA451	PA452	PA453	PA454	PA455	PA456	PA45AN
	CORE-VENT 5.3								
	MICRO-VENT 4.25								
	SCREW-VENT 4.7								
ASTRA TECH	4.5, 5.0, ST		AS501 (.5 mm)	AS502	AS503	AS504	AS505	AS506	AS5AN
LIFECORE BIOMEDICAL	RENOVA 3.75	SRS350	SRS351	SRS352	SRS353	SRS354	SRS355	SRS356	SRS35AN
NOBEL BIOCARE	SELECT, REPLACE 3.5								
	SELECT, REPLACE 4.3	SRS450	SRS451	SRS452	SRS453	SRS454	SRS455	SRS456	SRS45AN
	SELECT, REPLACE 5.0		SRS501	SRS502	SRS503	SRS504	SRS505	SRS506	SRS5AN
DENTSPLY FRIADENT	ANKYLOS 3.5, 4.5, 5.5		AK401 (1.65 mm)		AK402	AK403	AK404	AK405	AK4AN
	FRIALIT-2 & XIVE 3.4		FX341	FX342	FX343	FX344	FX345	FX346	FX34AN
	FRIALIT-2 & XIVE 3.8		FX381	FX382	FX383	FX384	FX385	FX386	FX38AN
	FRIALIT-2 & XIVE 4.5		FX451	FX452	FX453	FX454	FX455	FX456	FX45AN
	FRIALIT-2 & XIVE 5.5		FX551	FX552	FX553	FX554	FX555	FX556	FX55AN
LIFECORE BIOMEDICAL	RENOVA 4.5 ___		SW381	SW382	SW383	SW384	SW385	SW386	SW38AN
ZIMMER - SWISSPLUS	SWISSPLUS 3.8								
	SWISSPLUS 4.8		SW481	SW482	SW483	SW484	SW485	SW486	SW48AN
INNOVA	4.8 ENDOPORE (Internal Connection)								
LIFECORE BIOMEDICAL	STAGE-1								
ITI STRAUMANN	3.3, 4.1, 4.8								
BUCK MEDICAL RESEARCH	BMR								
OSSTEM	SSII 4.1 4.8								
OSSTEM	SSIII 4.0 5.0								
COWELL MEDI	ATLAS INTERNAL 3.5 4.0 5.0		T401	T401	T403	T404	T405	T406	T4AN
DENTIS	DENTIS INTERNAL 3.7 4.3 4.8								
MEGAGEN	EXFEEL INTERNAL 3.5 4.1 4.8								
NEOBIOTECH	SINUS QUICK IT 3.5 4.0 5.0								
	INTERNAL 3.3 4.1 4.8								
CHAORUM	INTERNAL 3.3 4.0 4.8								
DIO	INTERNAL PSI , IFI 3.5,4.1,4.8								
CAMLOG	CAMLOG 3.8		CAM381	CAM382	CAM383	CAM384	CAM385	CAM386	CAM38AN
	CAMLOG 4.3		CAM431	CAM432	CAM433	CAM434	CAM435	CAM436	CAM43AN
	CAMLOG 5.0		CAM501	CAM502	CAM503	CAM504	CAM505	CAM506	CAM5AN
	CAMLOG 6.0		CAM601	CAM602	CAM603	CAM604	CAM605	CAM606	CAM6AN
SYBRON IMPLANT SOLUTIONS	ENDOPORE INTERNAL HEX 4.1	IH410	IH411	IH412	IH413	IH414	IH415	IH416	IH41AN
	ENDOPORE INTERNAL HEX 5.0		IH511	IH512	IH513	IH514	IH515	IH516	IH51AN
	PITT-EASY 3.25/3.3		PE331	PE332	PE333	PE334	PE335	PE336	PE33AN
	PITT-EASY 3.75/3.8		PE381	PE382	PE383	PE384	PE385	PE386	PE38AN
	PITT-EASY 4.0/4.1		PE411	PE412	PE413	PE414	PE415	PE416	PE41AN
	PITT-EASY 4.9		PE491	PE492	PE493	PE494	PE495	PE496	PE49AN
IMPLANT DESCRIPTION	3.7 SCREWPLANT		SP371	SP372	SP373	SP374	SP375	SP376	SP37AN
	4.7 SCREWPLANT		SP471	SP472	SP473	SP474	SP475	SP476	SP47AN
	5.7 SCREWPLANT		SP571	SP572	SP573	SP574	SP575	SP576	SP57AN
IMTEC CORPORATION	4.3 ENDURE		ED431	ED432	ED433	ED434	ED435	ED436	ED43AN